



Immerse strip for 2 seconds. Remove with pads face up. SHAKE ONCE TO REMOVE EXCESS WATER. Read immediately: FCI>TCI>Alk>pH

Sumergir la tira durante 2 segundos. Sacarla del agua con las almohadillas hacia arriba. SACUDIR UNA VEZ PARA ELIMINAR EL EXCESO DE AGUA.

Comparer el color inmediatamente:
FCI>TCI>Alk>pH

Immerger la bande pour 2 secondes. Retirer la face au dessus. SECOUER POUR ÉLIMINER L'EXCÉDENT D'EAU. Lecture immédiate:
FCI>TCI>Alk>pH

LaMotte Company, Chestertown, MD 21620 USA www.lamotte.com/insta

Free Chlorine (Bromine) Cloro libre (bromo) Chlore libre (brome)

FCI (Br)	0	0.5	1	3	5	10
ppm	(0)	(1)	(2)	(6)	(10)	(20)

Total Chlorine Cloro total Chlore total

TCI	0	0.5	1	3	5	10
ppm						

Alkalinity Alcalinidad Alcalinité

Alk	0	40	80	120	180	240
ppm						

pH	6.2	6.8	7.2	7.8	8.4	9.0

Purple color indicates excessive sanitizer.
El tono morado indica un exceso de desinfectante.
Si violet, chlore ou brome trop élevé.

Code 2978

